


BLK Antibody / B lymphocyte kinase (R35515)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35515-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 640 |
| Applications | Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This BLK antibody is available for research use only. |



Western blot testing of human thymus lysate with BLK antibody at 0.1ug/ml. Predicted molecular weight: ~58 kDa, routinely observed at ~55 kDa.

Description

Additional name(s) for this target protein: B lymphoid tyrosine kinase; p55-Blk

Application Notes

Optimal dilution of the BLK antibody should be determined by the researcher.

Immunogen

Amino acids QWSPLKVSAQDKD were used as the immunogen for this BLK antibody.

Storage

Aliquot and store the BLK antibody at -20oC.