


BLNK Antibody / B-cell linker protein [clone 28B74] (RQ8578)

Catalog No.	Formulation	Size
RQ8578	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28B74
Purity	Affinity chromatography
UniProt	Q8WV28
Applications	Western Blot : 1:500-1:2000
Limitations	This BLNK antibody is available for research use only.



Western blot testing of human Raji cell lysate with BLNK antibody. Predicted molecular weight: ~51 kDa but can be observed at 65-70 kDa.

Description

Functions as a central linker protein, downstream of the B-cell receptor (BCR), bridging the SYK kinase to a multitude of signaling pathways and regulating biological outcomes of B-cell function and development. [UniProt]

Application Notes

Optimal dilution of the BLNK antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human BLNK protein was used as the immunogen for the BLNK antibody.

Storage

Store the BLNK antibody at -20oC.