


BLVRB Antibody (R36047)

Catalog No.	Formulation	Size
R36047-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	645
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This BLVRB antibody is available for research use only.



Western blot testing of human liver lysate with BLVRB antibody at 0.1ug/ml. Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: Biliverdin reductase B

Application Notes

Optimal dilution of the BLVRB antibody should be determined by the researcher.

Immunogen

Amino acids SRLPSEGPRPAH were used as the immunogen for this BLVRB antibody.

Storage

Aliquot and store the BLVRB antibody at -20oC.