


BMP15 Antibody (F54078)

Catalog No.	Formulation	Size
F54078-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54078-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O95972
Applications	Western Blot : 1:4000
Limitations	This BMP15 antibody is available for research use only.



Western blot testing of BMP15 antibody at 1:4000 + human testis lysate. Routinely observed molecular weight: ~50 kDa pro form and 16-17 kDa mature form.

Description

May be involved in follicular development. Oocyte- specific growth/differentiation factor that stimulates folliculogenesis and granulosa cell (GC) growth. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the BMP15 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Recombinant human BMP15 was used as the immunogen for the BMP15 antibody.

Storage

Aliquot the BMP15 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.