


BMP3 Antibody (F45244)

Catalog No.	Formulation	Size
F45244-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F45244-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P12645
Applications	Western Blot : 1:1000
Limitations	This BMP3 antibody is available for research use only.



Description

BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation.

Application Notes

Titration of the BMP3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 318-347 from the human protein was used as the immunogen for this BMP3 antibody.

Storage

Aliquot the BMP3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.