


BMP7 Antibody (R36451)

Catalog No.	Formulation	Size
R36451-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P18075
Gene ID	655
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This BMP7 antibody is available for research use only.



Western blot testing of pig kidney lysate with BMP7 antibody at 1ug/ml. Predicted molecular weight: ~49 kDa.

Description

Additional name(s) for this target protein: Bone morphogenetic protein 7

Application Notes

Optimal dilution of the BMP7 antibody should be determined by the researcher.

Immunogen

Amino acids 120-131 (SLQDSHFLTDAD) were used as the immunogen for this BMP7 antibody.

Storage

Aliquot and store the BMP7 antibody at -20oC.