


BMP7 Antibody (R36451)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R36451-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse, Rat, Dog, Pig, Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| UniProt | P18075 |
| Gene ID | 655 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000 |
| Limitations | This BMP7 antibody is available for research use only. |



Western blot testing of pig kidney lysate with BMP7 antibody at 1ug/ml. Predicted molecular weight: ~49 kDa.

Description

Additional name(s) for this target protein: Bone morphogenetic protein 7

Application Notes

Optimal dilution of the BMP7 antibody should be determined by the researcher.

Immunogen

Amino acids 120-131 (SLQDSHFLTDAD) were used as the immunogen for this BMP7 antibody.

Storage

Aliquot and store the BMP7 antibody at -20oC.