


BMPR2 Antibody (R33524)

Catalog No.	Formulation	Size
R33524-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	659
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BMPR2 antibody is available for research use only.



Western blot testing of human liver lysate with BMPR2 antibody at 1ug/ml. Predicted molecular weight ~115 kDa (precursor).

Description

Additional name(s) for this target protein: Bone morphogenetic protein receptor, type II

Application Notes

Optimal dilution of the BMPR2 antibody should be determined by the researcher.

Immunogen

Amino acids DLETNKLDPKEVDK were used as the immunogen for this BMPR2 antibody.

Storage

Aliquot and store the BMPR2 antibody at -20oC.