


# BMPR2 Antibody (R33524)

Catalog No.	Formulation	Size
R33524-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	659
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BMPR2 antibody is available for research use only.



#### **Description**

Additional name(s) for this target protein: Bone morphogenetic protein receptor, type II

## **Application Notes**

Optimal dilution of the BMPR2 antibody should be determined by the researcher.

### **Immunogen**

Amino acids DLETNKLDPKEVDK were used as the immunogen for this BMPR2 antibody. **Storage** Aliquot and store the BMPR2 antibody at -20oC.