


BPI Antibody (R34476)

Catalog No.	Formulation	Size
R34476-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	671
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This BPI antibody is available for research use only.



Western blot testing of human bone marrow lysate with BPI antibody at 0.3ug/ml.
Predicted molecular weight: ~54kDa.

Description

Additional name(s) for this target protein: Bactericidal permeability-increasing protein; BPIFD1, CAP 57

Application Notes

Optimal dilution of the BPI antibody should be determined by the researcher.

Immunogen

Amino acids KIPDYSDSFKIKH were used as the immunogen for this BPI antibody.

Storage

Aliquot and store the BPI antibody at -20oC.