


BPNT1 Antibody (R35703)

Catalog No.	Formulation	Size
R35703-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10380
Applications	Western Blot : 1-5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This BPNT1 antibody is available for research use only.



Western blot testing of human liver lysate with BPNT1 antibody at 1 ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: Bisphosphate 3'-nucleotidase

Application Notes

Optimal dilution of the BPNT1 antibody should be determined by the researcher.

Immunogen

Amino acids RNYDYYASRVESIK were used as the immunogen for this BPNT1 antibody.

Storage

Aliquot and store the BPNT1 antibody at -20oC.