


## Brain Creatine Kinase Antibody / CKB (R34577)

Catalog No.	Formulation	Size
R34577-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1152
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Brain Creatine Kinase antibody is available for research use only.



Western blot testing of human brain lysate with brain Creatine Kinase antibody at 0.1ug/ml. Predicted molecular weight: ~43 kDa.

### Description

Additional name(s) for this target protein: CKB; Creatine kinase B

### Application Notes

Optimal dilution of the Brain Creatine Kinase antibody should be determined by the researcher.

## **Immunogen**

Amino acids KPSDEHKTDLNPDN were used as the immunogen for this Brain Creatine Kinase antibody.

## **Storage**

Aliquot and store the Brain Creatine Kinase antibody at -20oC.