


## Brain Creatine Kinase Antibody / CKB (R34577)

Catalog No.	Formulation	Size
R34577-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Rat, Dog, Pig, Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	1152
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:2000
<b>Limitations</b>	This Brain Creatine Kinase antibody is available for research use only.



Western blot testing of human brain lysate with brain Creatine Kinase antibody at 0.1ug/ml. Predicted molecular weight: ~43 kDa.

### Description

Additional name(s) for this target protein: CKB; Creatine kinase B

This Creatine Kinase B Western Blot Antibody is part of a broader [Creatine Kinase B antibody for WB panel](#) offered by NSJ Bioreagents.

## **Application Notes**

Optimal dilution of the Brain Creatine Kinase antibody should be determined by the researcher.

## **Immunogen**

Amino acids KPSDEHKTDLNPDN were used as the immunogen for this Brain Creatine Kinase antibody.

## **Storage**

Aliquot and store the Brain Creatine Kinase antibody at -20oC.