


Bruno-like 5 Antibody (R33921)

Catalog No.	Formulation	Size
R33921-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60680
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Bruno-like 5 antibody is available for research use only.



Western blot testing of human heart lysate with Bruno-like 5 antibody at 1ug/ml.
Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: BRUNOL5; Bruno L5

Application Notes

Optimal dilution of the Bruno-like 5 antibody should be determined by the researcher.

Immunogen

Amino acids EPPGGQPDGMKD were used as the immunogen for this Bruno-like 5 antibody.

Storage

Aliquot and store the Bruno-like 5 antibody at -20oC.