

BST2 Antibody (R33878)

Catalog No.	Formulation	Size
R33878-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	684
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BST2 antibody is available for research use only.



Western blot testing of human spleen lysate with BST2 antibody at 0.3ug/ml. Expected molecular weight: 20-40 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Bone marrow stromal cell antigen 2, Tetherin, CD317

Application Notes

Optimal dilution of the BST2 antibody should be determined by the researcher.

Immunogen

Amino acids ELTEAQKGFQD were used as the immunogen for this BST2 antibody.

Storage

Aliquot and store the BST2 antibody at -20oC.