


BTBD7 Antibody (R33994)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33994-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Cow, Dog, Pig, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 55727 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This BTBD7 antibody is available for research use only. |



Western blot of MDCK lysate overexpressing mouse Btbd7 (lane 3) tested with BTBD7 antibody. Lane 1: mock transfection, lane 2: full construct chemically suppressed

Description

Additional name(s) for this target protein: BTB/POZ domain-containing protein 7

Application Notes

Optimal dilution of the BTBD7 antibody should be determined by the researcher.

Immunogen

Amino acids ESKLYSLDHGHEKPQ were used as the immunogen for this BTBD7 antibody.

Storage

Aliquot and store the BTBD7 antibody at -20oC.