


BTD Antibody (R33819)

Catalog No.	Formulation	Size
R33819-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	686
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This BTD antibody is available for research use only.



Western blot testing of human liver lysate with BTD antibody at 0.5ug/ml. Predicted molecular weight: ~61kDa.

Description

Additional name(s) for this target protein: Biotinidase

Application Notes

Optimal dilution of the BTD antibody should be determined by the researcher.

Immunogen

Amino acids EKDAQEVHCDEATK were used as the immunogen for this BTD antibody.

Storage

Aliquot and store the BTD antibody at -20°C.