


BTLA Antibody / CD272 (R34392)

Catalog No.	Formulation	Size
R34392-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q7Z6A9
Gene ID	151888
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This BTLA antibody is available for research use only.



Western blot testing of mouse spleen lysate with BTLA antibody at 1ug/ml. Predicted molecular weight ~33kDa.

Description

Additional name(s) for this target protein: B and T lymphocyte associated

Application Notes

Optimal dilution of the BTLA antibody should be determined by the researcher.

Immunogen

Amino acids KEAPTEYASICVRS were used as the immunogen for this BTLA antibody.

Storage

Aliquot and store the BTLA antibody at -20oC.