


c-Fos Antibody [clone RM374] (R20393)

Catalog No.	Formulation	Size
R20393-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM374
Purity	Protein A purified from animal origin-free supernatant
UniProt	P01100
Localization	Nuclear, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:500-1:1000 Western Blot : 1:1000-1:5000
Limitations	This recombinant c-Fos antibody is available for research use only.



Western blot testing of lysate from HeLa cells, non-treated (-) or treated (+) with TPA, with recombinant c-Fos antibody at 1:5000 dilution. Expected molecular weight: ~40 kDa (unmodified), 53-68 kDa (phosphorylated).



IHC staining of FFPE human thyroid tissue with recombinant c-Fos antibody at 1:1250.

Description

This antibody reacts to human Proto-oncogene c-Fos.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant c-Fos antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to N-terminus of human Proto-oncogene c-Fos was used as the immunogen for the recombinant c-Fos antibody.

Storage

Store the recombinant c-Fos antibody at -20oC.