


C14orf68 Antibody (R34651)

Catalog No.	Formulation	Size
R34651-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	283600
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This C14orf68 antibody is available for research use only.



IHC testing of FFPE human liver with C14orf68 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Chromosome 14 open reading frame 68

Application Notes

Optimal dilution of the C14orf68 antibody should be determined by the researcher.

Immunogen

Amino acids KSRLQADGQGQRR were used as the immunogen for this C14orf68 antibody.

Storage

Aliquot and store the C14orf68 antibody at -20oC.