

C6orf58 Antibody / LEG1 (F54531)


Catalog No.	Formulation	Size
F54531-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54531-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q6P5S2
Applications	Flow Cytometry : 1:25 (1x10 ⁶ cells) Immunohistochemistry (FFPE) : 1:25 Western Blot : 1:500-1:2000
Limitations	This C6orf58 antibody is available for research use only.




Western blot testing of human K562 cell lysate with C6orf58 antibody. Predicted molecular weight ~40 kDa.



Western blot testing of human liver tissue lysate with C6orf58 antibody. Predicted molecular weight ~40 kDa.



Flow cytometry testing of human K562 cells with C6orf58 antibody; Blue=isotype control, Green= C6orf58 antibody.



IHC testing of FFPE human testis carcinoma with C6orf58 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.

Description

C6orf58/LEG1 may be involved in early liver development. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the C6orf58 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 205-233 from the human protein was used as the immunogen for the C6orf58 antibody.

Storage

Aliquot the C6orf58 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.