


C7 Antibody / Complement C7 (R34789)

Catalog No.	Formulation	Size
R34789-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	730
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This C7 antibody is available for research use only.



Western blot testing of human liver lysate with C7 antibody at 0.1ug/ml. Predicted molecular weight: ~94 kDa.

Description

Additional name(s) for this target protein: Complement component 7

Application Notes

Optimal dilution of the C7 antibody should be determined by the researcher.

Immunogen

Amino acids NELTGQFRNRVINTK were used as the immunogen for this C7 antibody.

Storage

Aliquot and store the C7 antibody at -20oC.