


# C7 Antibody / Complement C7 (R34789)

Catalog No.	Formulation	Size
R34789-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	730
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This C7 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Complement component 7

### **Application Notes**

Optimal dilution of the C7 antibody should be determined by the researcher.

### **Immunogen**

Amino acids NELTGQFRNRVINTK were used as the immunogen for this C7 antibody.

# **Storage** Aliquot and store the C7 antibody at -20oC.