


C9 Antibody / Complement C9 (F43252)


Catalog No.	Formulation	Size
F43252-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43252-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P02748
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This C9 antibody is available for research use only.



C9 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human placenta tissue.



C9 antibody western blot analysis in ZR-75-1 lysate

Description

This gene encodes the final component of the complement system. It participates in the formation of the Membrane Attack Complex (MAC). The MAC assembles on bacterial membranes to form a pore, permitting disruption of bacterial membrane organization. Mutations in this gene cause component C9 deficiency. [provided by RefSeq].

Application Notes

Titration of the C9 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 191-220 from the human protein was used as the immunogen for this C9 antibody.

Storage

Aliquot the C9 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.