


C9 Antibody / Complement C9 (R34018)

Catalog No.	Formulation	Size
R34018-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	735
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This C9 antibody is available for research use only.



Western blot testing of human placenta lysate with C9 antibody at 2ug/ml. Predicted molecular weight: ~63/70kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Complement component C9

Application Notes

Optimal dilution of the C9 antibody should be determined by the researcher.

Immunogen

Amino acids ETKGEKNFRTEH were used as the immunogen for this C9 antibody.

Storage

Aliquot and store the C9 antibody at -20oC.