


CA1 Antibody (R35192)

Catalog No.	Formulation	Size
R35192-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	759
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CA1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Carbonic anhydrase I

Application Notes

Optimal dilution of the CA1 antibody should be determined by the researcher.

Immunogen

Amino acids QAIKTKGKRAP were used as the immunogen for this CA1 antibody. **Storage** Aliquot and store the CA1 antibody at -20oC.