


CA1 Antibody (R35192)

Catalog No.	Formulation	Size
R35192-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	759
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CA1 antibody is available for research use only.



Western blot testing of human liver lysate with CA1 antibody at 0.03ug/ml.
Predicted/observed molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: Carbonic anhydrase I

Application Notes

Optimal dilution of the CA1 antibody should be determined by the researcher.

Immunogen

Amino acids QAIKTKGKRAP were used as the immunogen for this CA1 antibody.

Storage

Aliquot and store the CA1 antibody at -20oC.