


CA12 Antibody Biotin Conjugate (R35918BTN)

Catalog No.	Formulation	Size
R35918BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	771
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:128000
Limitations	This CA12 antibody is available for research use only.



Western blot testing of human kidney lysate with biotinylated CA12 antibody at 0.1ug/ml.
Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Carbonic anhydrase XII, carbonic anhydrase 12

Application Notes

Optimal dilution of the CA12 antibody should be determined by the researcher.

Immunogen

Amino acids SHLQHVKYKGQE were used as the immunogen for this CA12 antibody.

Storage

Aliquot and store the CA12 antibody at -20oC.