


# CACNB2 Antibody (R33724)

Catalog No.	Formulation	Size
R33724-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	783
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This CACNB2 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Calcium channel, voltage-dependent, beta 2

## **Application Notes**

Optimal dilution of the CACNB2 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids ECNKQRSRHKSKDRY were used as the immunogen for this CACNB2 antibody. **Storage** Aliquot and store the CACNB2 antibody at -20oC.