


CACNB2 Antibody (R33724)

Catalog No.	Formulation	Size
R33724-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	783
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This CACNB2 antibody is available for research use only.



Western blot testing of Jurkat lysate with CACNB2 antibody at 1ug/ml. Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: Calcium channel, voltage-dependent, beta 2

Application Notes

Optimal dilution of the CACNB2 antibody should be determined by the researcher.

Immunogen

Amino acids ECNKQRSRHKSKDRY were used as the immunogen for this CACNB2 antibody.

Storage

Aliquot and store the CACNB2 antibody at -20oC.