


## Cadherin 23 Antibody (R36389)

Catalog No.	Formulation	Size
R36389-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64072
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Cadherin 23 antibody is available for research use only.



Western blot testing of human amygdala lysate with Cadherin 23 antibody at 1ug/ml.  
Predicted/observed molecular weight: ~60kDa.

## Description

Additional name(s) for this target protein: CDH23

## Application Notes

Optimal dilution of the Cadherin 23 antibody should be determined by the researcher.

## Immunogen

Amino acids YNISLYENVTVGTS were used as the immunogen for this Cadherin 23 antibody.

## **Storage**

Aliquot and store the Cadherin 23 antibody at -20oC.