


Calpain 1 Antibody (R34767)

Catalog No.	Formulation	Size
R34767-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	823
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Calpain 1 antibody is available for research use only.



Western blot testing of human kidney lysate with Calpain 1 antibody at 0.3ug/ml.
Predicted molecular weight: ~82kDa.

Description

Additional name(s) for this target protein: CAPN1

Application Notes

Optimal dilution of the Calpain 1 antibody should be determined by the researcher.

Immunogen

Amino acids NASRARSEQFINLR were used as the immunogen for this Calpain 1 antibody.

Storage

Aliquot and store the Calpain 1 antibody at -20oC.