


Calponin 3 Antibody (R34187)

Catalog No.	Formulation	Size
R34187-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1266
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Calponin 3 antibody is available for research use only.



Western blot testing of human duodenum lysate with Calponin 3 antibody at 0.3ug/ml.
Predicted molecular weight: ~36kDa.

Description

Additional name(s) for this target protein: Calponin, acidic; CNN3

Application Notes

Optimal dilution of the Calponin 3 antibody should be determined by the researcher.

Immunogen

Amino acids HGEYQDDYPRDYQYS were used as the immunogen for this Calponin 3 antibody.

Storage

Aliquot and store the Calponin 3 antibody at -20oC.