


CAPN1 Antibody (R35877)

Catalog No.	Formulation	Size
R35877-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	823
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CAPN1 antibody is available for research use only.



Western blot testing of human spleen lysate with CAPN1 antibody at 0.3ug/ml. Predicted molecular weight: ~82kDa.

Description

Additional name(s) for this target protein: Calpain 1

Application Notes

Optimal dilution of the CAPN1 antibody should be determined by the researcher.

Immunogen

Amino acids SEKSAGTVELDD were used as the immunogen for this CAPN1 antibody.

Storage

Aliquot and store the CAPN1 antibody at -20oC.