


Caspase 3 Antibody (R35575)

Catalog No.	Formulation	Size
R35575-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	836
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Caspase 3 antibody is available for research use only.



Western blot testing of MOLT4 lysate with Caspase 3 antibody at 1ug/ml. The expected ~32kDa band and the additional ~65kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: CASP3, SCA-1, YAMA, Apopain

Application Notes

Optimal dilution of the Caspase 3 antibody should be determined by the researcher.

Immunogen

Amino acids RDVSKEDHSKRS were used as the immunogen for this Caspase 3 antibody.

Storage

Aliquot and store the Caspase 3 antibody at -20oC.