


Caspase 4 Antibody (large subunit) (R33261)

Catalog No.	Formulation	Size
R33261-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cat, Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	837
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Caspase 4 antibody is available for research use only.



Western blot testing of human heart lysate with Caspase 4 antibody at 0.3ug/ml. Expected molecular weight ~43 kDa (precursor), ~30 kDa (pro domain + large subunit), ~20 kDa (large subunit).

Description

Additional name(s) for this target protein: CASP4; ICH-2

Application Notes

Optimal dilution of the Caspase 4 antibody should be determined by the researcher.

Immunogen

Amino acids CKERAEIYPIKERN were used as the immunogen for this Caspase 4 antibody.

Storage

Aliquot and store the Caspase 4 antibody at -20oC.