


Caspase 6 Antibody (R34902)

Catalog No.	Formulation	Size
R34902-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	839
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Caspase 6 antibody is available for research use only.



Western blot testing of human colon lysate with Caspase 6 antibody at 0.5ug/ml.
Predicted molecular weight: ~33 kDa (precursor).

Description

Additional name(s) for this target protein: Apoptotic protease MCH-2; CASP6

Application Notes

Optimal dilution of the Caspase 6 antibody should be determined by the researcher.

1. This Caspase 6 antibody is specific for isoform alpha.

Immunogen

Amino acids NRKVSQRRVDFCKDP were used as the immunogen for this Caspase 6 antibody.

Storage

Aliquot and store the Caspase 6 antibody at -20oC.