


Caspase 6 Antibody (R35479)

Catalog No.	Formulation	Size
R35479-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	839
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Caspase 6 antibody is available for research use only.



Western blot testing of human colon lysate with Caspase 6 antibody at 0.3ug/ml.
Predicted molecular weight: ~33 kDa (precursor).

Description

Additional name(s) for this target protein: CASP6, Apoptotic protease MCH-2

Application Notes

Optimal dilution of the Caspase 6 antibody should be determined by the researcher.

Immunogen

Amino acids QTEKLDTNITEVD were used as the immunogen for this Caspase 6 antibody.

Storage

Aliquot and store the Caspase 6 antibody at -20oC.