


Caspase-8 Antibody (small subunit) (R30599)


| Catalog No. | Formulation | Size |
|-------------|---|--------|
| R30599 | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Rabbit |
| Clonality | Polyclonal (rabbit origin) |
| Isotype | Rabbit IgG |
| Purity | Antigen affinity |
| Buffer | Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal |
| UniProt | Q14790 |
| Applications | Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml |
| Limitations | This Caspase-8 antibody is available for research use only. |



Western blot testing of Caspase-8 antibody and human HeLa cell lysate. Predicted molecular weight: ~55 kDa (pro), ~40 kDa (large + small subunit), ~11 kDa (small subunit).



IHC-P: Caspase-8 antibody testing of human breast cancer tissue

Description

The human CASP8 gene, whose product is also known as Caspase-8 and FLICE, encodes an interleukin-1beta converting enzyme (ICE)-related cysteine protease that is activated by the engagement of several different death receptors. It is immediately recruited to the Fas receptor once it oligomerizes, and its protease activity is crucial for the apoptotic response generated by the resulting death-inducing signaling complex (DISC). This gene contains at least 11 exons spanning approximately 30kb on human chromosome band 2q33-34. This region of chromosome 2 was previously reported as the location of the CASP10 gene, whose product is closely related to CASP8. Caspase-8 deficiency in humans is compatible with normal development and shows that the enzyme has a postnatal role in immune activation of naive lymphocytes.

Application Notes

The stated application concentrations are suggested starting points. Titration of the Caspase-8 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Amino acids 458-473 (TILTEVNYEVS NKDDK-human) were used as the immunogen for this Caspase-8 antibody. This sequence is from the small subunit.

Storage

After reconstitution, the Caspase-8 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.