


CBF Antibody (R33553)

Catalog No.	Formulation	Size
R33553-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12607
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CBF antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with CBF antibody at 0.5ug/ml. Predicted molecular weight: ~120kDa, observed here at ~140kDa.

Description

Additional name(s) for this target protein: CCAAT/enhancer binding protein zeta; Cebpz

Application Notes

Optimal dilution of the CBF antibody should be determined by the researcher.

Immunogen

Amino acids DNASF_KQLKWEAERD were used as the immunogen for this CBF antibody.

Storage

Aliquot and store the CBF antibody at -20°C.