


## CBL Antibody (R35586)

Catalog No.	Formulation	Size
R35586-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	867
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CBL antibody is available for research use only.



Western blot testing of NIH3T3 lysate with CBL antibody at 2ug/ml. Expected molecular weight: 100-120 kDa.

## Description

Additional name(s) for this target protein: Casitas B lineage lymphoma proto oncogene, RNF55, RING finger protein 55

## Application Notes

Optimal dilution of the CBL antibody should be determined by the researcher.

## **Immunogen**

Amino acids REFVSISSPAHVAT were used as the immunogen for this CBL antibody.

## **Storage**

Aliquot and store the CBL antibody at -20oC.