


## CBR3 Antibody (R35282)

Catalog No.	Formulation	Size
R35282-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	874
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CBR3 antibody is available for research use only.



Western blot testing of HEK293 lysate with CBR3 antibody at 0.5ug/ml. Predicted molecular weight: ~31kDa.

## Description

Additional name(s) for this target protein: Carbonyl reductase 3

## Application Notes

Optimal dilution of the CBR3 antibody should be determined by the researcher.

## Immunogen

Amino acids QGQLVHDKVVQNW were used as the immunogen for this CBR3 antibody.

## Storage

Aliquot and store the CBR3 antibody at -20oC.