


CCDC144A Antibody (F53974)

Catalog No.	Formulation	Size
F53974-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53974-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A2RUR9
Applications	Western Blot : 1:1000-2000
Limitations	This CCDC144A antibody is available for research use only.



Western blot testing of human 1) brain and 2) HL-60 cell lysate with CCDC144A antibody at 1:2000. Predicted molecular weight: 165/81/134 kDa (isoforms 1/2/3).

Description

Coiled coil domain containing 144A is a 1427 amino acid protein belonging to the CCDC family.

Application Notes

Titration of the CCDC144A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1257-1289 from the human protein was used as the immunogen for the CCDC144A antibody.

Storage

Aliquot the CCDC144A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.