


CCKAR Antibody (R34055)

Catalog No.	Formulation	Size
R34055-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	886
Applications	Immunohistochemistry (FFPE) : 10ug/ml ICC/IF : suitable ELISA (peptide) LOD : 1:4000
Limitations	This CCKAR antibody is available for research use only.



IHC testing of FFPE human pancreas with CCKAR antibody at 10ug/ml. HIER: microwaved with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Cholecystokinin A receptor

Application Notes

Optimal dilution of the CCKAR antibody should be determined by the researcher.

Immunogen

Amino acids FEASQKKSARKPST were used as the immunogen for this CCKAR antibody.

Storage

Aliquot and store the CCKAR antibody at -20oC.