


CCM2 Antibody (R34440)

Catalog No.	Formulation	Size
R34440-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	83605
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CCM2 antibody is available for research use only.



Western blot testing of human heart lysate with CCM2 antibody at 0.03ug/ml. Predicted molecular weight: ~49kDa.

Description

Additional name(s) for this target protein: Cerebral cavernous malformation 2

Application Notes

Optimal dilution of the CCM2 antibody should be determined by the researcher.

1. This CCM2 antibody is specific for isoforms 1 and 2/4.

Immunogen

Amino acids KGEKSRDKKAHEK were used as the immunogen for this CCM2 antibody.

Storage

Aliquot and store the CCM2 antibody at -20oC.