


CCNDBP1 Antibody (R36189)

Catalog No.	Formulation	Size
R36189-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23582
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This CCNDBP1 antibody is available for research use only.



Western blot testing of human placenta lysate with CCNDBP1 antibody at 0.3ug/ml.
Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Cyclin D-type binding-protein 1, DIP1, GCIP, MAID

Application Notes

Optimal dilution of the CCNDBP1 antibody should be determined by the researcher.

Immunogen

Amino acids TKASHVTPQPED were used as the immunogen for this CCNDBP1 antibody.

Storage

Aliquot and store the CCNDBP1 antibody at -20oC.