


## CCT3 Antibody (R34092)


Catalog No.	Formulation	Size
R34092-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7203
Localization	Cytoplasmic
Applications	Western Blot : 0.01-0.03ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CCT3 antibody is available for research use only.



Western blot testing of mouse testis lysate with CCT3 antibody at 0.01ug/ml. Predicted molecular weight: ~61 kDa.



IHC testing of FFPE human tonsil tissue with CCT3 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: Chaperonin containing TCP1, subunit 3

## Application Notes

Optimal dilution of the CCT3 antibody should be determined by the researcher.

## Immunogen

Amino acids GHKKKGDDQSRQGG were used as the immunogen for this CCT3 antibody.

## Storage

Aliquot and store the CCT3 antibody at -20oC.