


## CD20 Antibody [clone 407CT20.1.2] (F53598)

Catalog No.	Formulation	Size
F53598-0.1ML	In ascites with 0.09% sodium azide	0.1 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Ascites
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	407CT20.1.2
Purity	Ascites
UniProt	P11836
Applications	Western Blot : 1:1000-1:3200
Limitations	This CD20 antibody is available for research use only.



CD20 antibody western blot analysis in Ramos lysate.

### Description

CD20 may be involved in the regulation of B-cell activation and proliferation. [UniProt]

### Application Notes

Titration of the CD20 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 266-294 from the human protein was used as the immunogen for this CD20 antibody.

## **Storage**

Aliquot the CD20 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.