

CD200R1 Antibody / OX-2 receptor [clone CD200R1/9919R] (V5863)

Catalog No.	Formulation	Size
V5863-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5863-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5863SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	CD200R1/9919R
UniProt	Q8TD46
Localization	Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml Western Blot : 2-4ug/ml
Limitations	This CD200R1/OX-2 receptor antibody is available for research use only.

Application Notes

1. Optimal dilution of the CD200R1/OX-2 receptor antibody should be determined by the researcher.
2. This CD200R1/OX-2 receptor antibody is recombinantly produced by expression in CHO cells.

Immunogen

A recombinant fragment (around amino acids 50-250) of human CD200R1 protein (exact sequence is proprietary) was used as the immunogen for the CD200R1/OX-2 receptor antibody.

Storage

CD200R1/OX-2 receptor antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

