


CD22 Antibody / CD22 molecule [clone MSVA-022R] (V6145)

Catalog No.	Formulation	Size
V6145-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6145-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-022R
UniProt	P20273
Localization	Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This CD22/CD22 molecule antibody is available for research use only.



B-cell receptor CD22 Rabbit Recombinant Monoclonal Antibody (MSVA-022R) tested on many normal and cancer tissues. The immunohistochemistry staining in these tissues aligns with the expression data in Human Protein Atlas.

Manual Protocol: Freshly cut sections should be used (less than 10 days between cutting and staining). Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Target Retrieval Solution buffer. Apply the antibody at a dilution of 1:150 at 37°C for 60 minutes. Visualization of bound antibody by the EnVision Kit (Dako, Agilent) according to the manufacturer's directions.

Application Notes

- Optimal dilution of the CD22/CD22 molecule antibody should be determined by the researcher.
- This CD22/CD22 molecule antibody is recombinantly produced by expression in CHO cells.

Immunogen

A recombinant fragment (around amino acids 52-178) of human BL-CAM/CD22 protein (exact sequence is proprietary) was used as the immunogen for the CD22/CD22 molecule antibody.

Storage

CD22/CD22 molecule antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.